



## DK50 PLUS

Pre jednu stomatologickú súpravu  
Do jednego unitu stomatologicznego  
Для одной стоматологической установки



Príslušenstvo/ Akcesoria/ Принадлежности: **BOX MD KJF FS REG AOK**

	DK50 PLUS	DK50 PLUS/M	DK50 PLUS S	DK50 PLUS S/M
	1	1	1	1
DRYER	-	✓	-	✓
BOX	-	-	✓	✓
(5-7bar)	75 l.min <sup>-1</sup> / 5bar (50Hz) 85 l.min <sup>-1</sup> / 5bar (60Hz)	58 l.min <sup>-1</sup> / 5bar (50Hz) 68 l.min <sup>-1</sup> / 5bar (60Hz)	75 l.min <sup>-1</sup> / 5bar (50Hz) 85 l.min <sup>-1</sup> / 5bar (60Hz)	58 l.min <sup>-1</sup> / 5bar (50Hz) 68 l.min <sup>-1</sup> / 5bar (60Hz)
(6-8bar)	70 l.min <sup>-1</sup> / 6bar (50Hz) 80 l.min <sup>-1</sup> / 6bar (60Hz)	60 l.min <sup>-1</sup> / 6bar (50Hz) 70 l.min <sup>-1</sup> / 6bar (60Hz)	70 l.min <sup>-1</sup> / 6bar (50Hz) 80 l.min <sup>-1</sup> / 6bar (60Hz)	60 l.min <sup>-1</sup> / 6bar (50Hz) 70 l.min <sup>-1</sup> / 6bar (60Hz)
(8-10bar)	60 l.min <sup>-1</sup> / 8bar (50Hz) 70 l.min <sup>-1</sup> / 8bar (60Hz)	50 l.min <sup>-1</sup> / 8bar (50Hz) 60 l.min <sup>-1</sup> / 8bar (60Hz)	60 l.min <sup>-1</sup> / 8bar (50Hz) 70 l.min <sup>-1</sup> / 8bar (60Hz)	50 l.min <sup>-1</sup> / 8bar (50Hz) 60 l.min <sup>-1</sup> / 8bar (60Hz)
S1-S3	S1 - 100%	S1 - 100%	S1 - 100%	S1 - 100%
	≤ 66 dB(A)	≤ 66 dB(A)	≤ 47 dB(A)	≤ 47 dB(A)

\* Počet zubných súprav závisí od požiadaviek v jednotlivých krajinách.

\* Liczba unitów stomatologicznych zależy do potrzeb.

\* Количество стоматологических установок зависит от региональных требований.

## DK50 2V

Pre dve stomatologické súpravy  
Do dwóch unitów stomatologicznych  
Для двух стоматологических установок



Príslušenstvo/ Akcesoria/ Принадлежности: **BOX MD KJF FS REG AOK**

	DK50 2V	DK50 2V/M	DK50 2VS	DK50 2VS/M
	2	2	2	2
DRYER	-	✓	-	✓
BOX	-	-	✓	✓
(5-7bar)	140 l.min <sup>-1</sup> / 5bar	108 l.min <sup>-1</sup> / 5bar	140 l.min <sup>-1</sup> / 5bar	108 l.min <sup>-1</sup> / 5bar
(6-8bar)	135 l.min <sup>-1</sup> / 6bar	104 l.min <sup>-1</sup> / 6bar	135 l.min <sup>-1</sup> / 6bar	104 l.min <sup>-1</sup> / 6bar
(8-10bar)	110 l.min <sup>-1</sup> / 8bar	87 l.min <sup>-1</sup> / 8bar	110 l.min <sup>-1</sup> / 8bar	87 l.min <sup>-1</sup> / 8bar
S1-S3	S1 - 100%	S1 - 100%	S1 - 100%	S1 - 100%
	≤ 67 dB(A)	≤ 69 dB(A)	≤ 51,5 dB(A)	≤ 55 dB(A)

\* Počet zubných súprav závisí od požiadaviek v jednotlivých krajinách.

\* Liczba unitów stomatologicznych zależy do potrzeb.

\* Количество стоматологических установок зависит от региональных требований.